





NIPRO EXTENSION TUBE

DÂY NỐI TRUYỀN DỊCH/ DÂY NỐI BƠM TIÊM ĐIỆN



- Dây nối làm bằng PVC y tế, Non-Toxic, Non-Pyrogenic
- Nhiều lựa chọn chiều dài khác nhau
- Lựa chọn đầu nối dạng thẳng, có nắp cố định
- Lựa chọn đầu nối vặn xoắn chống rò rỉ, dây có khóa áp lực dương chống trào ngược dịch
- Tiệt trùng ETO, đạt chuẩn ISO 13485

NIPRO EXTENSION TUBE

Hình sản phẩm	Mã sản phẩm	Mô tả	Đóng gói
	C-NN-WET-3010-ECE	- Dây dài 300mm - Đầu nối dạng thẳng, có nắp cố định	50pcs/hộp
	C-NN-WET-5025-ECE	- Dây dài 450mm - Đầu nối dạng thẳng, có nắp cố định	50pcs/hộp
	C-NN-WET-9950-ECE	- Dây dài 950mm - Đầu nối dạng thẳng, có nắp cố định	50pcs/hộp
	C-NN-WET-1500-E	- Dây dài 1500mm - Đầu nối vặn xoắn chống rò rỉ, dây có khóa áp lực dương chống trào ngược dịch	100pcs/hộp

*Dùng 1 lần
Hạn dùng 5 năm
ISO13485*

NIPRO CORPORATION

3-9-9, Honjo-Nishi, Kita-ku, Osaka,
531-8510 Japan
Phone: +81 6 6372 2331
Fax: +81 6 6371 7422

NIPRO SALES (THAILAND) CO., LTD.

89/170, Nipro Building, Vipavadee Rangsit Road,
Talat Bangkhen, Laksi, Bangkok 10210, Thailand
Phone: +66 2551 2000 Fax: +66 2551 2006
Nst@niprosales.com

NIPRO SALES VIETNAM REPRESENTATIVE OFFICE

38A Le Van Huan, Tan Binh District,
HCMC, Vietnam
Phone: +84 932702380
Nhuy@niprosales.com

Manufacturer: NIPRO (THAILAND) CORPORATION LIMITED

10/2 Moo 8 Bangnomko, Sena Phra Nakorn Si Ayutthaya, Thailand 13110

www.nipro.co.jp/en/index.html

www.niprosales.com