




# NIPRO STOPCOCK


An toàn-Dễ kết nối-Hiệu suất chịu áp lực cao



- Chịu áp lực cao, kháng lipid và tương thích với nhiều loại dịch truyền
- Trong suốt dễ dàng quan sát dòng chảy
- Dễ dàng kết nối chắc chắn với các dụng cụ khác
- Sản xuất tại Nhật Bản
- Tiệt trùng: khí E.O

# NIPRO STOPCOCK



MODEL	ORDER CODE	PICTURE	MATERIAL	SPECIFICATION
3W-R - L	A-NCN-3W-R-L-E		<b>Body: Polycarbonate</b> <b>Lever: Polyethylene</b>	2 đầu nối female luer lock vặn xoắn 1 đầu nối male luer lock vặn xoắn Khóa xoay 360 độ

Hạn dùng: 5 năm  
 Sử dụng 1 lần  
 ISO13485, EC

*Các thiết kế của Nipro ra đời từ vô số ý tưởng và công nghệ đổi mới.  
 Chất lượng và hiệu suất vượt trội từ nhiều năm kinh nghiệm kiểm soát sản xuất nghiêm ngặt.*

**NIPRO CORPORATION**  
 3-9-9, Honjo-Nishi, Kita-ku, Osaka,  
 531-8510 Japan  
 Phone: +81 6 6372 2331  
 Fax: +81 6 6371 7422

**NIPRO SALES (THAILAND) CO., LTD.**  
 89/170, Nipro Building, Viphavadee Rangsit Road,  
 Talad Bangkhen, Laksi, Bangkok 10210, Thailand  
 Phone: +66 2551 2000 Fax: +66 2551 2006  
 Nst@niprosales.com

**NIPRO SALES VIETNAM  
 REPRESENTATIVE OFFICE**  
 38A Le Van Huan, Tan Binh District,  
 HCMC, Vietnam  
 Phone: +84 932702380  
 Nhuy@niprosales.com

Manufacturer

**NIPRO CORPORATION ODATE FACTORY**

8-7, Hanukiyachi, Niida, Odate-shi, Akita, 018-5794 Japan  
[www.nipro.co.jp/en/index.html](http://www.nipro.co.jp/en/index.html)  
[www.niprosales.com](http://www.niprosales.com)