



NIPRO SPINAL NEEDLES

KIM CHỌC DÒ/ GÂY TÊ TỦY SỐNG NIPRO

Chính xác- An toàn- Tin cậy



- Công nghệ thép sản xuất tại Nhật Bản, thiết kế chính xác giúp định hướng đúng mặt vát và dễ dàng thâm nhập vào vị trí khó
- Thành kim mỏng cho dòng chảy tốt
- Có khoang trống giữa que thông nòng và thành kim, chuỗi kim có các rãnh ngang giúp cầm chắc chắn, trong suốt giúp phát hiện nhanh dịch não tủy chảy ra
- Mã code theo màu sắc để chọn lựa kích thước kim
- Tiệt trùng: EO


NIPRO SPINAL NEEDLE

*Các thiết bị y tế của Nipro ra đời từ vô số ý tưởng và công nghệ đổi mới.
Chất lượng và hiệu suất vượt trội từ nhiều năm kinh nghiệm và kiểm soát sản xuất nghiêm ngặt.*




Bevel geometry

Lancet Needle



Pencil Needle



Nhiều lựa chọn

Code	Gauge(mm)	Color Code	Needle Length
E-SPN-T1690E	16G (1.65mm)	White	1-1/2"(38mm) or 2-1/2"(60mm) or 3"(75mm) or 3-1/2"(88mm)
E-SPN-T1890E	18G (1.25mm)	Pink	
E-SPN-T1990E	19G (1.05mm)	Cream	
E-SPN-T2090E	20G (0.90mm)	Yellow	
E-SPN-T2190E	21G (0.80mm)	Deep Green	
E-SPN-T2290E	22G (0.70mm)	Black	
E-SPN-T2390E	23G (0.63mm)	Deep Blue	
E-SPN-T2590E	25G (0.50mm)	Orange	
E-SPN-T2690E	26G (0.46mm)	Brown	
E-SPN-T2790E	27G (0.41mm)	Grey	

*Hạn dùng 5 years
Sử dụng một lần
ISO13485, EC*

NIPRO CORPORATION

3-9-9, Honjo-Nishi, Kita-ku, Osaka,
531-8510 Japan
Phone: +81 6 6372 2331
Fax: +81 6 6371 7422

NIPRO SALES (THAILAND) CO., LTD.

89/170, Nipro Building, Viphavadee Rangsit Road,
Talat Bangkhen, Laksi, Bangkok 10210, Thailand
Phone: +66 2551 2000 Fax: +66 2551 2006
Nst@niprosales.com

NIPRO SALES VIETNAM REPRESENTATIVE OFFICE

38A Le Van Huan, Tan Binh District,
HCMC, Vietnam
Phone: +84 932702380
Nhuy@niprosales.com

Manufacturer

NIPRO MEDICAL INDUSTRIES LTD. TATEBAYASHI PLANT

2-19-64, Matsubara, Tatebayashi-shi, Gunma, 374-8518 Japan
www.nipro.co.jp/en/index.html
www.niprosales.com